

INSTRUCTIONS FOR USE

DENTAL DESENSITIZER

1. PRODUCT DESCRIPTION

Dental Desensitizer AF do BRASIL is a low viscosity desensitizing gel based on potassium nitrate, and has a neural desensitizing action generated by Potassium Nitrate. Potassium nitrate has been considered a very effective agent in the treatment of hyperesthetic dentin. Its action is based on the fact that the increase in extracellular potassium concentration depolarizes nerve fiber membranes, blocking the passage of stimuli, reducing hypersensitivity. Must be used for professional use only.

2. PRESENTATION

1. Dental Desensitizer 2.0%

- In polypropylene syringes containing 2g, 2.5g, 3g, 4g, 5g or 10g.
- In single packs or containing 2 units; 3 units; 4 units; or 10 pcs.

2. Dental Desensitizer 0.2%

- In polypropylene syringes containing 2g, 2.5g, 3g, 4g, 5g or 10g.
- In single packs or containing 2 units; 3 units; 4 units; or 10 pcs.

3. COMPOSITION

- Active Ingredients: Potassium Nitrate and Sodium Fluoride.
- Inactive Ingredients: Deionized Water, Glycerin, neutralizing and thickening agent.

4. INDICATION OF USE

Indicated in the treatment of teeth with low, moderate and high sensitivity and was developed to relieve dental sensitivity triggered by thermal and chemical agents.

It is also indicated in the treatment of sensitivity that can be caused by tooth whitening.

5. PRODUCT ACTION

It is indicated for use before, during and after tooth whitening and when sensitivity is significant.

6. CONSERVATION/ STORAGE



STORE IN A DRY PLACE AT ROOM TEMPERATURE(15 TO 30°C).

7. INSTRUCTIONS FOR USE

- Clean the surface of the teeth by removing plaque, calculus and other deposits that may prevent the product from penetrating the teeth;
- Apply the desensitizing gel evenly over the teeth and wait for a minimum of 10 minutes;
- Alternatively, the gel can be applied using the bleaching tray, applying it as if it were the bleach itself;
- After the application time, remove the gel from the teeth with cotton and plenty of water. Prevent the patient from swallowing the desensitizer;
- Repeat the application of the desensitizer as needed or in the case of tooth whitening, following the treatment time. The number of applications should not exceed once a day.

Note: After the desensitizing time has elapsed and with it still on the teeth, it is suggested to rub the product with a micro-applicator, brush or felt discs at low speed and intermittent circular movements (low pressure).

This procedure mechanically assists the penetration of the desensitizing agent into the tooth, thus improving its performance.

8. CONTRAINDICATIONS

Do not use on patients with any sensitivity to the components of this product.

9. WARNINGS

- Avoid intentional and excessive ingestion of the product;
- If any problem occurs during or after application, notify your dentist;
- Non-injectable product;
- Do not reuse packaging.

10. SPECIAL CARES

If a large volume of product is accidentally swallowed, seek medical attention.

11. MANIPULATION

Aseptic environment.

12. INTERACTIONS

There is not.