

INSTRUCTIONS FOR USE

PREFORMED SECTIONAL STAINLESS STEEL MATRIX ASSORTED KIT

1. PRODUCT DESCRIPTION

The Preformed Sectional Stainless Steel Matrix Kit AF do BRASIL for use in restorations, its ease and agility due to being pre-cut reduces operative time, having a good adaptation and providing more comfort to the patient.

2. PRESENTATION

Assorted Kit of preformed sectional stainless steel dies:

- 02 Fixing Clamps;
- 08 Silicone Protectors;
- 80 Assorted Matrices S, M, L and XL.

3. COMPOSITION

Stainless steel and Silicone.

4. INDICATION OF USE

Indicated for dental use to aid in class II restorations.

5. PRODUCT ACTION

Preformed sectional stainless steel matrix allows a better natural anatomy of the tooth after the procedure.

6. CONSERVATION / STORAGE

- Keep the product in its packaging always closed;
- Store the Product in a cool dry place at a temperature below 60°C. Humidity should be less than 95° RH

7. INSTRUCTIONS FOR USE

When starting the procedure, check the size of the preformed sectional stainless steel matrix (S, M, L and XL) that best adapts to the cavity to be restored.

Place the matrix on the proximal surface to be restored, then the wedge and then the fixation clip with the silicone protectors.

Note: The use or not of the clamp and silicone protectors and the way they are installed depend on the user's personal preferences and the clinical situation. Perform the restore process. Once completed, remove the clamp, then the wedge, and finally the die.

8. CONTRAINDICATIONS

There is not.

9. WARNINGS

- Do not ingest;
- Keep out of the reach of children;
- For Dental use only;
- Handle according to the instructions;
- Sterilize before use.

10. SPECIAL CARES

Professional use, single use, non-sterile, discard after use.

11. MANIPULATION

aseptic environment

12. INTERACTIONS

There is not.