

INSTRUCTIONS FOR USE

BUCAL MIRROR

Clinical Odontological Mirror

1. PRODUCT DESCRIPTION

The AF do BRASIL Clinical Mirror is an autoclavable product, manufactured in stainless steel, which brings high quality, reliability and resistance. Intended for dental professionals who seek efficiency and quality in the care of their patients.

2. PRESENTATION

Box containing 12 units and/or individual;

Available in the following models:

- Clinical Mirror Plan nº 3;
- Clinical Mirror Plan nº 5;
- Clinical Mirror 1st Plan Front Surface nº 5.

3. COMPOSITION

Stainless Steel 304 and Glass.

4. INDICATION FOR USE

It is indicated as an image reflector, used in the visualization of the dental arch, in interventions and support in dental procedures.

5. STORAGE

Store the product in a cool, dry place, free from direct sunlight.

6. INSTRUCTIONS FOR USE

- The mirror must be attached (screwed) to the handle that is sold separately;
- Take it into the oral cavity so that it can be used in the procedures intended.

7. CONTRA INDICATIONS

Not applicable.

8. WARNINGS

Always check that the product is in perfect condition and free of defects;

9. SPECIAL CARE

Not applicable.

10. VALIDITY

AF do BRASIL Clinical Mirror is valid for 5 years after the manufacturing date.