



Efficiency, Quality and Fair Value.



## Light-Curing Gingival Barrier

The Gengival Barrier AF do BRASIL is a light-curing resin that does not irritate the patient's gingival tissues. It has adequate viscosity and thixotropy to allow easy and perfect coverage of soft tissues, without flowing into undesirable regions. It has excellent sealing power and minimal and perfect adhesion on the gum. It comes out completely when you take it out. Yields up to 5 full arcade applications.

## INDICATION

Indicated as a substitute for a rubber sheet to protect the gingival tissue and in cases of risk of irritation from contact with products used in clinical procedures, such as tooth whitening. It also acts as an adjunct in absolute isolation for better sealing of a cavity preparation or endodontic treatment. Protect the product from direct sunlight and store at temperatures between 5 and 28°C.

## **PRESENTATION**

Box with 1 (one) syringe containing 2g of Light-Curing Gengival Barrier AF do BRASIL + 2 (two) applicator tips.



nº 80432410063







