



**AF<sup>do</sup>BRASIL**  
Produtos Odontológicos

Efficiency, Quality and Fair Value.

# Dental Radiological Fixer and Developer



The AF do BRASIL Dental Radiological Fixer is indicated for image fixation on the film. The Fixer, together with the Developer, allows the image extracted from the patient to be transferred to the film so that the exam can be viewed in any place that has light. Ensures faster and longer-lasting development.

## INDICATION

Indicated for fixing the developed image on dental x-ray film.

## PRESENTATION

Bottle containing 500ml of Fixer AF do BRAZIL.

Bottle containing 500ml of Developer AF Dental do BRAZIL.



📍 Rua Damaris Rosa de Menezes Monteiro, 136  
Pq. Tec. Francisco Sciarra | Lindóia | 86031-216 | Londrina/PR  
📞 @afdobrasil 📧 @afdobrasil



+55 43 3158.2886  
sac@afdobrasil.ind.br  
www.afdobrasil.ind.br