



AF do **BRASIL**
Produtos Odontológicos

Efficiency, Quality and Fair Value.

Dental Radiological Fixer and Developer



The Dental Radiologic Developer of AF do BRASIL is used to develop the exam image on the film film. The Developer together with the Fixer allows the image extracted from the patient to be transferred to the film so that the exam can be viewed in any place that has light. Ensures faster development with greater durability. **READY TO USE.**

INDICATION

Indicated for fixing the developed image on dental x-ray film.

PRESENTATION

Bottle containing 500ml of Fixer AF do BRASIL.

Bottle containing 500ml of Developer AF Dental do BRASIL.



📍 Rua Damaris Rosa de Menezes Monteiro, 136
Pq. Tec. Francisco Sciarra | Lindóia | 86031-216 | Londrina/PR
📞 @afdobrasil 📧 @afdobrasil



+55 43 3158.2886
sac@afdobrasil.ind.br
www.afdobrasil.ind.br